

StatLIA® Immunoassay Troubleshooting Chart

Match the pattern of (Flagged Parameters) for each Component on the chart with the Flagged Parameters (p<0.01) on the StatLIA Report to pinpoint the failed Component. Flags with a cross-hatched pattern are only listed with those parameters used to gauge whether a problem with a Component has significantly affected the results. These parameters may or may not be flagged, depending on the severity of the problem. (Parameters along the x-axis appear in the same order as they do in the "Qutlier" tables at the top of the StatLIA Report.)



The Following Parameters are Used to Gauge Whether a Potential Problem With One of the Components has Significantly Affected the Results



^{*} Monitor Points include: Minimum/Maximum Detectable Concentrations, ED20, ED80, %Bmax/Tracer Activity, %NSB/Tracer Activity. Curve Probability includes: Parallel Probability and Compare Probability. (All of these parameters are computed automatically by StatLIA.)